


Disposable OPU Needle for veterinary use  
cow ova vacuuming needle ...

# COVA Needle

Registered No.715 Ministry of Agriculture, Forestry and Fisheries of Japan



Needle Point  
(Blast treated)

## Needle Set for Ultrasound Guided Ovum Pickup for Cow

- Applicable Ultrasound Probes : Aloka, Honda, Hitachi and Others
- Individually Packed
- Sterilized by Ethylene Oxide Gas
- 25 pcs/box




## Our COVA Needles are...

- One-way type
- Single use only : Free from washing, drying, sterilizing...
- Aspiration tube : Already assembled 1500mm length PVC tube.
- Needle size : 17G x 500mm (ALOKA), 17G x 490mm (HONDA)  
17G x 540mm (HITACHI), 17G x 600mm (C-Type)

## What's more...

- Blast treated needle point : For clear monitor image.
- Already Assembled, Packaged and Sterilized : Just open and use it.
- All parts (Needle, Handle, and Tube) are bonded firmly to ensure the prevention of separation.
- Needles are coated by silicone oil (same as medical injection needles) for smooth and delicate handling.

see the charts



## Caution :

- Single use only - disposable product
- Do not use if the pack is opened or damaged.
- For safety and sanitation, discard the product after use.
- Store in cool, dry, dark place.

MISAWA MEDICAL INDUSTRY CO., LTD.

351, Asahi-machi Kasama-shi, Ibaraki-ken, 309-1717 Japan

Phone : +81-296-77-8804 , Fax : +81-296-77-8849

email : sales@misawa-medical.co.jp http://www.misawa-medical.com